


Aksesoret e CT-Skaner model “Ingenuity Flex 32” ose ekuivalent

- 1.mbjtëse e kokës së sheshtë
- 2.Shiringe Injeksioni (Shiringes)te kontrastit
- 3.Mbjtësin e Injeksionit(Shiringes)te Kontrastit
- 4.Paneli i kontrollit te Shiringes
- 5.Rrypat e sigurimit te Pacientit
- 6.Printer

Pjeset qe permban Skaneri

- 1.Gantry
- 2.Shtratin e pacientit
- 3.UPS
- 4.Isotran LM(Transformatorin)



Minimum room layout
285 sq.ft. (26.5 sq. meters)

This layout represents minimum space for equipment installation and service access.
National and local codes may require additional space for electrical safety and patient access.

Dimensions and weights*

	Weight	Height	Width	Depth
1 Gantry	1950 kg (4300 lbs.)	203 cm (79")	239 cm (94")	94 cm (37")
2 Table	404 kg (890 lbs.)	107 cm (42")	69 cm (27")	515 cm (203")
3 Operator console	88 kg (193 lbs.)	118 cm (47")	120 cm (47")	91 cm (36")
4 Computer/server (CRC)	142 kg (313 lbs.)	77 cm (31")	61 cm (24")	89 cm (35")
5 UPS for CRC	136 kg (300 lbs.)	46 cm (18")	63 cm (25")	66 cm (26")